RSID FIELD KITS

Designed for Crime Scene Investigation

Designed to withstand extreme temperatures (-23°C to 48°C)

The RSID™- Saliva, Semen and Blood Field Kits are designed for fast, easy and reliable detection of specific human fluids from a variety of samples encountered at crime scenes.

Kits don't require refrigeration or special temperature storage. As these are used for crime scene investigations, they are designed to withstand extreme temperatures (-23°C to 48°C).

Components of the Kit:

- 10 Strip Tests RSID™Saliva or RSID™Semen or RSID™Blood
- 10 Extraction/Running buffer tubes pre-loaded
- 10 Transfer Pipettes Disposable single-use
- Cotton swabs
- DNA Free Water for wetting swabs
- Stain sampling/cutting protocol and forensic sample extraction method

RSID™ SALIVA FIELD KIT - 10 Tests /Kit Cat.° 1000-10 RSID™ SEMEN FIELD KIT - 10 Tests /Kit Cat.° 2000-10 RSID™ BLOOD FIELD KIT - 10 Tests /Kit Cat.° 3000-10









Tests

Multifluid Field Kit

SALIVA, SEMEN, BLOOD 15 tests/kit (5 of each). Each est comes with with pre-loaded buffer, swab, and

Multi Fluid FIELD KIT:

RSID™ Multi Fluid FIELD KIT - 15 Tests/kit Cat. 5500

Include:

Distributed by:

- 5 RSID-Saliva + 5 RSID-Semen + 5 RSID-Blood
- Universal Buffer pre-loaded
- DNA Free Water for wetting swabs
- Cotton swabs
- Transfer Pipettes

OCIETÀ ITALIANA CHIMICI Life Sciences

RSID" SID" SID" Via Rio nell'Elba, 140 - Rome / ITALY

Phone +39 06-8818936 / 8800211 - Fax +39 06-8815319 Email info@sichim.com Web-Site www.sichim.com